

## Beginning:

- First, we need to prepare the anatomical data. Since we have 4 T1w images, we need to preprocess all T1w images and then generate an averaged anatomical T1w image according to those 4 images.
- To do this, two methods are possible: 1) using ANT and 2) using SPM-12
- All preprocessing methods in here are for glasslexical task but the same process would be applied for motor task which consists of two runs.

## ANT averaging method:

```
1 import os
2 import ants
3
4 # Define paths to the four T1-weighted images
5 # in this sample code it will only happens for the subject number one, you can run that for each of your subject
6 t1_images = [
7     "D:\\sussex\\ids pro\\research group\\DATA\\raw data\\data\\ds000224_R1.0.2\\sub-MSC01\\ses-struct01\\anat\\sub-MSC01_ses-struct01_run-01_T1w.nii",
8     "D:\\sussex\\ids pro\\research group\\DATA\\raw data\\data\\ds000224_R1.0.2\\sub-MSC01\\ses-struct01\\anat\\sub-MSC01_ses-struct01_run-02_T1w.nii",
9     "D:\\sussex\\ids pro\\research group\\DATA\\raw data\\data\\ds000224_R1.0.2\\sub-MSC01\\ses-struct02\\anat\\sub-MSC01_ses-struct02_run-01_T1w.nii",
10    "D:\\sussex\\ids pro\\research group\\DATA\\raw data\\data\\ds000224_R1.0.2\\sub-MSC01\\ses-struct02\\anat\\sub-MSC01_ses-struct02_run-02_T1w.nii"
11 ]
12
13 # Load the first image as reference
14 ref_img = ants.image_read(t1_images[0])
15
16 # Initialize an empty list to store registered images
17 registered_imgs = [ref_img]
18
19 # Register all images to the reference
20 for img_path in t1_images[1:]:
21     img = ants.image_read(img_path)
22     # Supervisor should decide between "Rigid", "Affine", and "SyN" but "Affine" seems the best
23     reg_img = ants.registration(fixed=ref_img, moving=img, type_of_transform='Affine')['warpedmovout']
24     registered_imgs.append(reg_img)
25
26 # Compute the average image
27 avg_t1 = sum(registered_imgs) / len(registered_imgs)
28
29 # Save the averaged image
30 output_path = "D:\\sussex\\ids pro\\research group\\DATA\\T1w_avg.nii.gz"
31 ants.image_write(avg_t1, output_path)
32
33 print(f"Averaged T1-weighted image saved at: {output_path}")
```


## SPM-12 averaging method:

- four anatomical images need to be segmented to generate deformation fields. Then, using anatomical images and their deformation fields, these four images will be wrapped into MNI space.
- Then, 4 anatomical images will be smoothed using a 6 mm Gaussian kernel.
- The averaged anatomical image will be generated by averaging 4 normalized anatomical images (w prefixes) using the ImCalc option.

## Anatomical segmentation:


- The following image demonstrates the segmentation process. Images that need to be segmented should be added (e.g., in our case four anatomical images including 2 sessions and 2 runs in each session).
- The deformation field at the end of this batch folder should be set as "forward".

- This will result in a file including “deformation fields” for all 4 images. These images have “y\_” prefixes.




## Anatomical normalization:

- Normalization batch will be used in this part for normalizing all 4 anatomical images in the same MNI space.
- Since the deformation field was set to forward in the previous step, only “write” will be selected for this step of normalization and not “estimate and write”.
- For each anatomical image, the image will be added in “images to write” and its corresponding deformation field (y\_) will be added to the “deformation field” section.
- To do this for all 4 images, “replicate: subject” will be utilized while the same process is executing for the other three images.
- Files with “w” prefixes will be generated after normalization.




## Anatomical smoothing:

- Finally, the smoothing batch will be used to generate smoothed images from 4 anatomical images.
- All images will be added together.
- Files with “s” prefixes will be generated after smoothing.




## Anatomical averaging:


- ImClac batch folder will be selected, this batch will be used for averaging a number of images.




- In the input images, 4 normalized images will be selected including images with “w” prefixes from 2 sessions and their 2 runs (not adding images with “sw” prefixes, however, deciding the use of which image depends on your research-related factors and decision, in some cases “sw” images will be used).




- The output file's name will be selected (e.g., T1w\_averaged\_sub01)



- The output directory will be set in a path that you prefer.




- In the expression section, you need to add the formula for your averaging. As is presented in the following image, since 4 images exist, the sum of all images was divided by the number of images.



## CONN preprocessing pipeline:


# 90 subjects, including 10 sessions for each 9 subjects (each subject's session would be treated as a single subject).

# Since we have two runs for each glasslexical and motor, 2 sessions would be selected for each subject.




# Averaged structural images for each subject would be imported.


# the averaged structural image is the same for every 10 images (e.g., 1-10,11-20,...)



# The first and second runs of the glasslexical functional image will be imported to the first and second sessions. (The same thing happens for the motor task.)







```
# all masks will be imported in ROI section
```







# the existed condition will be deleted. The new condition will be imported as following:




# Which approach to follow? Having a whole condition named glasslexical or break this condition into its sections (e.g., noun, verb, and etc)




- The preprocessing will be selected from the toolbar located down the screen.




- According to the JSON file, interleaved (Simens) will be selected as slice order.



- All other parts would be set as default (only the Gaussian kernel will be set to 6mm)



- After preprocessing, the denoising part starts by clicking on Done.



- The denoising process will be shown. By clicking on “Done”, your data will be denoised and ready for the first-level analysis.

